

LabCollector

"Your Laboratory management solution"

User's guide

I-COLLECTOR HUB

Version 4.0 – February 2022



Summary

1) Installation of I-Collector Hub	1
1.1) <i>Configuration of the IP address of I-Collector Hub.....</i>	<i>1</i>
1.2) <i>Configuration of a Proxy server in I-Collector Hub.....</i>	<i>3</i>
1.3) <i>Checking the connection to I-Collector Server or LabCollector....</i>	<i>5</i>
1.4) <i>Connection of an instrument to I-Collector Hub for Instruments</i>	<i>6</i>
1.5) <i>Connection of AgileBio Data Logger to I-Collector Hub for Data Logger</i>	<i>9</i>
1.6) <i>Configuration of I-Collector Hub with a local I-Collector Server with an Internet connection.....</i>	<i>12</i>
1.7) <i>Configuration of I-Collector Hub with a local I-Collector Server without Internet connection</i>	<i>12</i>
4) Usage	13
5) Assistance	17
6) Specifications	20
6.1) <i>Hardware.....</i>	<i>20</i>
6.2) <i>Data buffering</i>	<i>21</i>
6.3) <i>Software Updates.....</i>	<i>22</i>
6.4) <i>Connection protocols.....</i>	<i>22</i>
6.5) <i>Connection limits.....</i>	<i>23</i>
7) Support	23
8) Version.....	23

1) Installation of I-Collector Hub

To install I-Collector Hub follow these steps:

- Configure the IP address of I-Collector Hub
- Configure a Proxy Server in I-Collector Hub if necessary
- Configure I-Collector Hub to connect to I-Collector Server or LabCollector
- Configure the instruments or data logger to connect to I-Collector Hub

1.1) Configuration of the IP address of I-Collector Hub

Important

I-Collector Hub can only work with a static IP address that must be defined into the DHCP server of your network.

I-Collector Hub won't work with a dynamic IP address assigned by a local DHCP server that changes with time.

The MAC address of each I-Collector Hub (like 9E-B6-D0-F2-85-73 for example) is indicated on the sticker on the back of I-Collector Hub. You can use this MAC address to associate one fixed IP address in your DHCP server to the corresponding I-Collector Hub.

1) To configure I-Collector Hub IP address:

- In your local network identify a free and fixed IP address to be assigned to I-Collector Hub. This address must not change in the future.
- Use the MAC address of one I-Collector Hub to associate a static IP address in the DHCP server to this Hub.

-
- 2) Plug I-Collector Hub to an Ethernet cable connected to your local network where data from instruments or loggers are accessible. I-Collector Server (local or in the Cloud) should be accessible from the local network where you connect I-Collector Hub.

Important:

I-Collector Hub should:

- have a direct access (or through a Proxy server) to Internet if I-Collector server or LabCollector is hosted in the Cloud,
- If the connection to Internet goes through a Proxy server and I-Collector server or LabCollector is hosted in the Cloud, start I-Collector Hub and go to paragraph “Configuration of a Proxy server in I-Collector Hub” to configure I-Collector Hub proxy,
- I-Collector Hub should always have access to I-Collector server or LabCollector if it is hosted in your local network.


- 3) Plug I-Collector to the power adaptor to start it.

If I-Collector Hub is connected to I-Collector Server through the local network and/or Internet and/or Proxy, and is already configured, it should automatically download the appropriate configuration from I-Collector Server to connect to instruments.

1.2) Configuration of a Proxy server in I-Collector Hub


If I-Collector Hub must use a proxy server access to access Internet, configure the proxy server in I-Collector Hub:

- Configure the IP address of the Hub (see paragraph “Configuration of the IP address of I-Collector Hub”)
- Wait for the IP address to be assigned and note it (WWW.XXX.YYY.ZZZ),
- In a web browser (Chrome, Edge, Firefox) connect to the IP address of the I-Collector Hub at <http://WWW.XXX.YYY.ZZZ:5000> to access the I-Collector Hub configuration page,




- Enter the username: hub
- Enter the password: THE_MAC_ADRESS_OF_THE_HUB (example: A0-CE-C8-39-8C-2C)

-
- Click Login to enter
 - In the proxy section, enter the proxy settings of the proxy server, click Save



The image shows a web form titled "Proxy Settings (Optional)". It contains three input fields: "Proxy URL", "Proxy Username", and "Proxy Password". Each field has a placeholder text with the same name as the field. Below the fields is a blue "Save" button.

- In the Server Status section click Stop and Start (to restart the internal services of the Hub and apply the new configuration)




The image shows a web interface section titled "Server Status". It displays the MAC Address as "B8-27-EB-01-7B-D0". Below this is a green progress bar labeled "Running". At the bottom, there are two buttons: a red "Stop" button and a green "Start" button.

I-Collector Hub should now connect to the Internet.

1.3) Checking the connection to I-Collector Server or LabCollector

To check if the connection between I-Collector Hub and I-Collector Server is up and running, connect to the I-Collector Hub administration page:

- In a web browser (Chrome, Edge, Firefox) connect to the IP address of the Hub at <http://WWW.XXX.YYY.ZZZ:5000> to access the I-Collector Hub configuration page
- Enter the username: hub
- Enter the password: THE_MAC_ADRESS_OF_THE_HUB (example: A0-CE-C8-39-8C-2C)
- Click Login to enter
- In the Connections List section you should see the connections established by the I-Collector Hub



- In the Server Logs section you should see the Logs of all the transactions going through I-Collector Server

Server Logs

Show entries Search:

	Name	Datetime	Message
1	system	2/16/2022 4:49...	0 files to transfer...
2	system	2/16/2022 4:49...	0 messages to transfer...
3	system	2/16/2022 4:49...	0 files to transfer...
4	system	2/16/2022 4:49...	0 messages to transfer...
5	system	2/16/2022 4:48...	0 files to transfer...
6	system	2/16/2022 4:48...	0 messages to transfer...
7	system	2/16/2022 4:48...	0 files to transfer...
8	system	2/16/2022 4:48...	0 messages to transfer...
9	system	2/16/2022 4:48...	0 files to transfer...
10	system	2/16/2022 4:48...	0 messages to transfer...

Showing 1 to 10 of 100 entries

First Previous 2 3 4 5 ... 10 Next Last

1.4) Connection of an instrument to I-Collector Hub for Instruments

Instruments with TCP/IP connection

If the instrument is directly connected to I-Collector Hub through a TCP/IP connection in the same network, indicate in the instrument administration panel the IP address (WWW.XXX.YYY.ZZZ) of the I-Collector Hub and the Port of the connection corresponding to the instrument Name found in the Connections List section of the administration page of I-Collector Hub.

Connections List							
	Name	Communication	Port	SFTP Port	SFTP Username	SFTP Password	Target Server
1	DataLogger	TCP Server	80				http://dev1.agilebio.com/neos...

Instrument with SFTP features

If the instrument can send result files on an SFTP server, the result files should be sent to I-Collector Hub as an SFTP server with the IP address of the I-Collector Hub (WWW.XXX.YYY.ZZZ), the Port of the connection corresponding to the instrument Name found in the Connections List section of the administration page of I-Collector Hub, the SFTP Username and Password indicated.

Connections List							
	Name	Communication	Port	SFTP Port	SFTP Username	SFTP Password	Target Server
1	DataLogger	TCP Server	80				http://dev1.agilebio.com/neos...


Instrument without SFTP features

If the instrument can only send result files in the local network, you must use and configure the I-Collector Client App to watch the folder and send the result files to I-Collector Hub.

To download I-Collector Client App connect to the LabCollector client area on <https://labcollector.com/>

To configure I-Collector Client App please read <https://labcollector.com/support/knowledge-base/how-do-i-install-and-configure-i-collector-client-app/>


I-Collector Client App must be configured to send the result files produced by the instrument to the I-Collector Hub acting as an SFTP server using the correct IP address of the I-Collector Hub (WWW.XXX.YYY.ZZZ) and the Port of the connection corresponding to the instrument Name found in the Connections List section of the administration page of I-Collector Hub.



	Name	Communication	Port	SFTP Port	SFTP Username	SFTP Password	Target Server
1	DataLogger	TCP Server	80				http://dev1.agilebio.com/neos...

Instruments with USB connection


If the instrument can only be connected by USB port, you should use an adapter to convert the USB to TCP/IP connection. Use the configuration page of the USB to TCP/IP adapter itself to send the result produced by the instrument to the I-Collector Hub via direct TCP/IP connection to the correct IP address of the I-Collector Hub (WWW.XXX.YYY.ZZZ) and the Port of the connection corresponding to the instrument Name found in the Connections List section of the administration page of I-Collector Hub.



	Name	Communication	Port	SFTP Port	SFTP Username	SFTP Password	Target Server
1	DataLogger	TCP Server	80				http://dev1.agilebio.com/neos...

Instruments with Serial connection

If the instrument can only be connected by Serial (RS232) port, you should use an adapter to convert the Serial(RS232) to TCP/IP connection. Use the configuration page of the Serial to TCP/IP adapter itself to send the result produced by the instrument to the I-Collector Hub via direct TCP/IP connection, using the correct IP address of the I-Collector Hub (WWW.XXX.YYY.ZZZ) and the Port of the connection corresponding to the instrument Name found in the Connections List section of the administration page of I-Collector Hub.




	Name	Communication	Port	SFTP Port	SFTP Username	SFTP Password	Target Server
1	DataLogger	TCP Server	80				http://dev1.agilebio.com/neos...

1.5) Connection of AgileBio Data Logger to I-Collector Hub for Data Logger

To configure the connection of an AgileBio Data Logger to a I-Collector Hub:

- In a web browser (Chrome, Edge, Firefox) enter the IP address of the Data Logger to connect to the configuration panel of the Data Logger
- Go to the menu System Configuration -> POST Client



- In the URL field indicate the IP address of I-Collector Hub without any port
- o Important: Uncheck POST Header

POST Client Configuration

Use this page to configure the automatic sending of data to an HTTP server using POST.

Enable Client:
 POST Header:
 Data Type: XML
 URL: http://192.168.0.241
 Period: 30 seconds
 Enable Proxy:
 Proxy URL: http://myproxy.com:8080
 Filter: Off (send all data)
 Device: All devices
 ROMId: 0000000000000000
 Offset: 0
 Tag:
 Tag:
 Tag:
 Tag:
 Tag:
 Tag:
 Tag:
 Tag:
 Tag:
 Tag:
 No Attributes:

- Click Save

In the administration page of I-Collector Hub you should see only one connection in the Connection List

Connections List							
	Name	Communication	Port	SFTP Port	SFTP Username	SFTP Password	Target Server
1	DataLo...	TCP Server	80				http://dev...


1.6) Configuration of I-Collector Hub with a local I-Collector Server with an Internet connection

If the local I-Collector Server has been installed by AgileBio, I-Collector Hub will automatically download from the Cloud the correct configuration to connect to the local I-Collector Server.

1.7) Configuration of I-Collector Hub with a local I-Collector Server without Internet connection

In the Connections List section of the administration page of I-Collector Hub:

- Indicate the Local Target URL corresponding to the right location of I-Collector Server on your network
- Click Save



The screenshot displays the 'Connections List' section of the I-Collector Hub administration interface. It features a table with the following data:

	Name	Communication	Port	SFTP Port	SFTP Username	SFTP Password	Target Server
1	DataLogger	TCP Server	80				http://dev1.agilebio.com/neos...

Below the table, the 'Local Target URL' field is highlighted with a green oval. The URL entered is `http://c.../UserFiles/lclogger_autoload_ow.php`. A red warning message below the field reads: 'Please do not change this value unless you are sure about it'. A blue 'Save' button is located at the bottom of the configuration area.


Note: If the local I-Collector Server has been installed by AgileBio, the provided I-Collector Hub may be shipped already configured.

4) Usage


I-Collector Hub is used to connect instruments or data loggers locally and send data to a LabCollector instance or I-Collector server (locally or in the Cloud) acting as a security HTTP/HTTPS bridge and buffering data;

I-Collector Hub is made for uni-directional communication only:

- I-Collector Hub is sending data from Instruments to I-Collector Server
- I-Collector Hub can't send data from I-Collector Server to the instrument because generally I-Collector Server can't connect directly to I-Collector Hub through a firewall.




I-Collector Hub is sending HTTPS data




I-Collector Hub is can store data locally in case the connection to I-Collector Server is not working


There are two types of I-Collector Hub:

- I-Collector Hub for Instruments
- I-Collector Hub for Data logger





I-Collector Hub for Data Logger



Example of usage of I-Collector Hub with I-Collector server in the Cloud

The main benefits of using I-Collector Hub are:

-
- Full Automation; when correctly configured I-Collector Hub can:

- connect to I-Collector server via Internet
- retrieve the connection configurations of your instruments from I-Collector server based on client ID defined in I-Collector server.
- Auto start receiving connection and data from instruments

- Buffering data

- I-Collector hub receive data (from data logger or instrument) and send them to LabCollector/I-Collector Server
- When Internet is not available it stores data (for a week maximum before deletion) and send them to LabCollector/I-Collector Server when Internet is back
- The buffering function secures the connection between the source of data (instruments, logger) and the target (I-Collector Server/LabCollector)
- I-Collector hub stores buffered data in a very powerful and stable database server (MySQL)
- I-Collector hub avoid data loss from instruments/loggers due to bad Internet connection

- Security

- I-Collector Hub receive data locally with HTTP protocol
- I-Collector Hub send data (to I-Collector Server / LabCollector) with secure HTTPS
- I-Collector Hub acts as a SECURITY BRIDGE for instruments to Internet

-
- I-Collector hub stores buffered data in a very powerful and stable database server (MySQL)
 - Ease of use of I-Collector Hub on Raspberry PI
 - I-Collector Hub on Raspberry PI is small
 - I-Collector Hub on Raspberry PI needs very little maintenance
 - Raspberry PI hardware is easy to replace because I-Collector Hub is contained on an SD Card
 - I-Collector Hub on Raspberry PI is easy to update because it is forward/backward compatible with Raspberry PI (from v3 <-> v4)

5) Assistance

I-Collector Hub doesn't start

In case I-Collector Hub does not start (no light), unplug and plug the power supply again. If I-Collector Hub does not start, check the power supply:


- If you have a spare I-Collector Hub, try the spare power supply.
- If the I-Collector Hub does not start with the spare power supply, contact our support

Checking the connection of I-Collector Hub to I-Collector Server / LabCollector

To check if I-Collector Hub is sending data to I-Collector Server/LabCollector, connect to the I-Collector Hub administration panel and check the Server Logs section:

- The list of the last 100 messages sent to I-Collector Server is indicated.

- The Date of the message is indicated in the Datetime column



The screenshot shows a window titled "Server Logs". At the top, there is a "Show 10 entries" dropdown and a "Search:" text box. Below this is a table with the following columns: "Name", "Datetime", and "Message". The table contains 10 rows of log entries, all from the "system" source. The "Datetime" column shows a sequence of timestamps from 2/16/2022 4:49:00 to 2/16/2022 4:48:00. The "Message" column contains alternating messages: "0 files to transfer..." and "0 messages to transfer...". Below the table, it says "Showing 1 to 10 of 100 entries" and a pagination bar with buttons for "First", "Previous", "1", "2", "3", "4", "5", "...", "10", "Next", and "Last".


	Name	Datetime	Message
1	system	2/16/2022 4:49...	0 files to transfer...
2	system	2/16/2022 4:49...	0 messages to transfer...
3	system	2/16/2022 4:49...	0 files to transfer...
4	system	2/16/2022 4:49...	0 messages to transfer...
5	system	2/16/2022 4:48...	0 files to transfer...
6	system	2/16/2022 4:48...	0 messages to transfer...
7	system	2/16/2022 4:48...	0 files to transfer...
8	system	2/16/2022 4:48...	0 messages to transfer...
9	system	2/16/2022 4:48...	0 files to transfer...
10	system	2/16/2022 4:48...	0 messages to transfer...

I-Collector Hub does not send data to I-Collector Server / LabCollector

If I-Collector Hub does not send data to I-Collector or LabCollector:

- Check I-Collector Hub is plugged and started,
- Check I-Collector Hub is connected to your local network by sending a ping request to its IP address for example,
- Check the instruments and/or data logger(s) are sending data to the network data repository (directory) or I-Collector Hub,


- Check I-Collector Hub is connected to instruments or data logger(s),
- Check I-Collector Hub is correctly configured to receive data from instruments or data logger(s),
- Check I-Collector Hub can connect and send data to I-Collector Server, through your local network or through Internet.
- In your LabCollector installation, restart the Plugin corresponding to the instrument:
 - o Connector to the LabCollector Parser add-on
 - o Go to the menu File Parsing History -> I-Collector viewer



- o Select the instrument that is not correctly connected



- Click the submenu Connection -> Restart Plugin



- This will restart the Plugin of I-Collector Server dedicated to the connection of the instrument

6) Specifications

6.1) Hardware

I-Collector Hub is installed on a tiny Linux computer called Raspberry PI <https://www.raspberrypi.com/>:

- *Raspberry PI 3* (1 Gb memory) for I-Collector Hub for Data Logger
- *Raspberry PI 4* (4 Gb memory) for I-Collector Hub for Instruments

Raspberry PI contains:

- 1 or more HDMI port to connect a display
- 2 or more USB ports to connect a keyboard, or can be used for update if network is unavailable
- 1 RJ45 Ethernet port to connect to the network
- 1 Wifi receiver/emitter to connect to the local Wifi
- 1 power port (220/110 V) (USB-C or micro-USB)

-
- 1 SD-Card of 16 GB as a disk

6.2) Data buffering

I-Collector Hub is buffering data received from instruments/logger in a MySQL database. In MySQL, data are deleted by I-Collector Hub one week after reception (data rotation).

In case of correct connection to I-Collector Server/LabCollector :

- Data are temporary stored in the MYSQL database
- Data are sent by batch every 10s to I-Collector Server/LabCollector

In case of interrupted connection to I-Collector Server/LabCollector:

- I-Collector Hub waits for the return of the connection to send data
- Then all pending data are sent to I-Collector Server/LabCollector

Maximum volume of data buffer

The size of the database buffering data is limited by the size of the disk of I-Collector Hub (SD Card in case of Raspberry PI).

Example of usage:

- For 2 dataloggers connected to 1 I-Collector Hub
- The 2 dataloggers sending one record every 5 seconds,
- The dataloggers sending data for 1 week
- => The MySQL database will be 500 MB after one week.

6.3) Software Updates

There is no planned update for I-Collector Hub on Raspberry PI. Updates will be performed by hardware replacement if necessary.

6.4) Connection protocols

HTTP/HTTPS connection

Instruments and data loggers connect by HTTP to I-Collector Hub. I-Collector Hub is connected by HTTPS to I-Collector Server / LabCollector

Instrument with Direct TCP/IP connection

Connection of the instrument to I-Collector Hub is done by direct TCP/IP connection.

Instrument with Serial (RS232) connection

Connection of the instrument to I-Collector Hub is done by direct TCP/IP connection with an RS232->TCP/IP adapter.

Instrument with USB connection

Connection of the instrument to I-Collector Hub is done by direct TCP/IP connection with an USB->TCP/IP adapter.

Result File on the network

Connection of the instrument to I-Collector Hub is done by using I-Collector Client App (<https://labcollector.com/support/knowledge-base/how-do-i-install-and-configure-i-collector-client-app/>) to send data from the network to I-Collector HUB acting as an SFTP server.

6.5) Connection limits

I-Collector Hub for Instrument

One I-Collector Hub for Instrument can be connected to a maximum of five instruments.

One I-Collector Hub for Instrument can't be connected to a Data Logger.

I-Collector Hub for Instrument v1 can act as an SFTP server but not as an FTP server to receive files.

I-Collector Hub for Data Logger

One I-Collector Hub for Data Logger can be connected to an unlimited number of Data Loggers.

One I-Collector Hub for Data Logger can't be connected to an instrument.

7) Support

To contact support please connect to LabCollector.com and open a ticket in your client area.

8) Version

I-Collector Hub user manual	
Version of the manual:	4.0
Date:	21 February 2022
Manual applies to version:	I-Collector Hub v1

Product name:	I-Collector Hub for Data Loggers
AgileBio product reference:	ICOL-HUB-LOG
Product name:	I-Collector Hub for Instruments
AgileBio product reference:	ICOL-HUB-INS



<http://www.labcollector.com>
sales@agilebio.com

AgileBio USA
5473 Kearny Villa Road Suite 255
San Diego, CA 92123
USA
Tel: 347 368 1315
Fax: (800) 453 9128

<http://www.agilebio.com>

AgileBio Headquarters
75 rue de Lourmel
75015 Paris
FRANCE
Tel: 01 41 79 15 85
Fax: 01 72 70 40 22